



24-Hour Urine Collection Container

Storing and measuring urine specimen during 24 hours



P.I.P. Urine Containers are made from polyethylene and are graduated to provide fast and easy measurement. These containers have two apertures: one is wide in order to facilitate sample collection by patient; the other is small for laboratory personnel to pour the specimen

into other containers without trickling. They are supplied in 2 models: horizontal and upright, with different capacities in two colors. The upright model occupies small space and the horizontal model is suitable for placing in refrigerators.

- Graduated
- Handle for easy transportation
- Separate sampling and transferring apertures
- Amber tinted model for photosensitive analytes



GS1 Code	Capacity	Model	Graduation (mm)	Opening Diameter	Color	Dimensions (mm)	Qty./pack
6260807500330	2.5 ml	Horizontal	250	φ 75	Colorless	240 × 115 × 160	30
6260807501061	3.5 ml	Horizontal	250	φ 75	Colorless	240 × 115 × 200	36
6260807501559	2.5 ml	Horizontal	250	φ 75	Amber tinted	240 × 115 × 160	30
6260807501566	3.5 ml	Horizontal	250	φ 75	Amber tinted	240 × 115 × 200	36
6260807501498	2.5 ml	Upright	100	φ 75	Colorless	130 × 130 × 230	48
6260807501504	3.5 ml	Upright	100	φ 75	Colorless	130 × 130 × 290	40
6260807501535	2.5 ml	Upright	100	φ 75	Amber tinted	130 × 130 × 230	48
6260807501542	3.5 ml	Upright	100	φ 75	Amber tinted	130 × 130 × 290	40

* For more information on the chemical resistance of polystyrene (PS), refer to our website at www.medpip.com.



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

Postal code:1576635714

P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159

Tel: 88765847 - 88765915 - 88545922-9

info@medpip.com

www.medpip.com

@poleideal

@poleideal

+98912 334 0197

Pole Ideal Pars Co.