



Pole Ideal Pars

(25-well)

Conical tube tray



Technical Information

These racks are designed in order to store and organize 15ml and 50ml conical tubes. One of the advantages of these racks is that a multitude of them can be connected together. They are also made of polypropylene, autoclavable, freezable, and resistant to most of the common laboratory acids and chemicals.

In order to increase efficiency, organization and differentiation of racks, they are available in 4 different colors (yellow, red, blue, and green).

Capacity	Conical Tube	Dimensions* (Length x width x Height)	Qty/pack
25-well	15ml	145 x 110 x 35	10
25-well	50ml	210 x 180 x 35	10

* millimeter



Pole Ideal Pars

Conical Tube



Technical Information

This disposable tube is made of transparent polypropylene, autoclavable, freezable, and resistant to most of the common laboratory acids and chemicals. The graduated conical tubes are designed in two sizes of 15 and 50 ml in order to store samples and also be used in centrifuge. They have screw cap and conical bases. These tubes can also be provided in Conical tube tray (25-well). A multitude of these racks can be connected together. In order to easily Classifying samples, racks and tube caps are available in 4 different colors (yellow, red, blue, and green).

Capacity	Graduation	Height *	Aperture External Diameter	Qty/pack
15ml	0.5ml	120	Ø 17	25
50ml	2.5ml	120	Ø 30	25

* millimeter