



Pole Ideal Pars Co.
www.medpip.com

Diamond Scribe

Marking on glass and hard surfaces, slides in particular

Diamond Scribe is for marking on hard surfaces such as slide, glass, etc. Its body is made from stainless steel and is available in two colors of silver and navy blue.

GS1 Code	Length
6260807500415	135 mm



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

Postal code:1576635714

P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159

Tel: 88765847 - 88765915 - 88545922-9

info@medpip.com

www.medpip.com

@poleideal

@poleideal

+98912 334 0197

Pole Ideal Pars Co.