



Disposable Base Mould

Paraffin embedding and preparing tissue block

P.I.P. Disposable Base Mould is designed with 6 different models with excellent thermal exchange and the paraffin does not stick to it, so that the expensive and time-consuming process of cleaning metal base moulds will be no longer required. The curved edges and the completely flat and smooth inner surface

allow easy block removal and paraffin ribbon continuity. Their built up sidewalls reduce risk of paraffin spillage during embedding and the standard dimensions are compatible with most of the tissue cassettes. They are cost-effective to be discarded after use, yet strong enough to be reused.

GS1 Code	Model	Inner well dimensions	Outer dimensions
6260807502013	PIP-07	7 mm × 7 mm × 6 mm	50 mm × 37 mm × 11 mm
6260807501993	PIP-15	15 mm × 15 mm × 6 mm	50 mm × 37 mm × 11 mm
6260807501979	PIP-24	24 mm × 24 mm × 6 mm	50 mm × 37 mm × 11 mm
6260807502006	PIP-30	30 mm × 24 mm × 6 mm	50 mm × 37 mm × 11 mm
6260807502020	PIP-37	37 mm × 24 mm × 6 mm	50 mm × 37 mm × 11 mm
6260807502242	PIP-30M	30 mm × 20 mm × 12 mm	55 mm × 39 mm × 18 mm

