



Filing Cabinets for Slide

Managing and archiving 76 mm × 26 mm slides

Because of the probable need to reuse a slide or preparation of new slides in hospitals pathology departments and different laboratory sections, long-term preserving and archiving is of great importance. With regard to the necessity and sensitivity of retrieving the samples, P.I.P. Filing Cabinets have the following properties:

- Managing and archiving 76 mm × 26 mm slides
- Humidity resistant in comparison to metal models
- No decay or mould in comparison to wooden models
- Lightweight in comparison to metal models
- High-capacity drawers with numbering for easy retrieval
- Smoothly-moving drawers
- No sharp edges
- Stackable
- Made from ABS



GS1 Code	Capacity	Outer dimensions (mm)
6260807502303	1200 slides	510 × 458 × 130

