



Pole Ideal Pars Co.  
www.medpip.com

## First Aid Kit

### Providing first aid requirements

Regarding the growing number of accidents that lead to injury and burns, it is even more necessary to have first aid kits available in any place. P.I.P. First Aid Kit can be used at home, office, school, workshop, gym, inside vehicles, and many other places. It is made from semitransparent polypropylene which is available in three sizes with different colors. As for the size and capacity, each model can be used to meet a wide range of requirements

and with different internal arrangements for specified applications such as burns, cuts, poisoning, asphyxia, sting, bleeding, hypothermia, eye care etc. In order to facilitate classifying instruments inside the kit, there are internal dividers for medium and large models. The lid has two latches in order for preventing unwanted opening of the lid. There is also a molded-in handle for easy transport.



GS1 Code	Model	Dimensions (mm)	Qty./Pack
6260807500378	S (Small)	95 × 30 × 125	6
6260807500385	M (Medium)	206 × 49 × 205	2
6260807500392	L (Large)	270 × 102 × 290	1



1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran  
Postal code:1576635714 P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159  
Tel: 88765847 - 88765915 - 88545922-9 info@medpip.com www.medpip.com  
@poleideal @poleideal ☎ +98912 334 0197 Pole Ideal Pars Co.