



## Hematology Cell Counter

### Manual differential counter



P.I.T. Hematology Cell Counter is made according to the international standards. It has an ergonomic keyboard which consists of 10 keys as Poly, Mono, Eos, nRbc, Baso, Band, Meta, Myeto, L.V.F, and Lymph for required parameters in Diff function. This device has audio alarm (after counting 100 white blood cells) and the percentage of required

parameters can be calculated simply by pressing the related key. In addition, counting stops automatically after reaching 100 cells. Therefore, in the event of pressing a key mistakenly, it will not disrupt the calculations. Yet, user can get the reticulocyte count and the actual number of leukocytes using the secondary keys %True Retic and True WBC.

Technical Data	Display	LCD- 16 × 2 character
	Display dimensions	62 mm × 13 mm
	Voltage	5 V
	Weight	365 gr

GS1 Code	Cat No.	Dimensions
6265738000010	RC902	205 mm x 135 mm x 55 mm

Unit11, No.13, Hojjatdoust st., South Felestin St., Keshavarz Blvd., Tehran, Iran

Tel:+98 21 88994482-3 Fax:+98 21 88994484 Postal Code:1416694334 Sales Office: +98 21 88545922