



Pole Ideal Pars

ROLLER MIXER



P.I.P. Hematology Mixer

P.I.P. Hematology Mixer is designed and produced according to the international standards to be used in medical diagnosis laboratories, research and medical centers

Scope of use of Roller Mixer:

- Mixing the blood containing anticoagulant and keeping it homogenize in order to prevent blood cells from becoming lysis and deformation
- Dissolving the lyophilized vials such as serum, controls, and calibrators vials
- Homogenizing and keeping suspensions (such as paint, oil, cosmetics and hygiene) used in research centers
- Preventing the sedimentation of the blood cells that is being tested

Technical Information:

Power: 220V, 50Hz
Speed: 60rpm

Dimensions: 50cm x 26cm x 12cm
Rollers Length: 32cm



FOR SAFETY BOX

1st floor, No.12, Naghdi St., Jahantab St., Mottahari Ave. Tehran, Iran
Tel: +98 21 88545922-9 / Fax: +98 21 88765561 & +98 21 88767159
Postal code: 1576635714 P.O. Box: 15875-9483