



Pole Ideal Pars Co.  
www.medpip.com


## Inoculating Loop Holder

Holding metal inoculating loops for sampling

To extend the application of Inoculating Loop Holder and its lifetime, it is strile, made from brass alloy which is resistant to heat. The handle of this holder is insulated by PVC

coating which increases user's convenience and prevents it from heating up. Each pack contains 6 holders plus 12 meters of special wire.

GS1 Code	Volume	Dimensions	Qty./pack
6260807500316	20 g	178 mm × $\phi$ 5 mm	6



1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.