



Pole Ideal Pars Co.  
www.medpip.com

## Integral Wash Bottle

Transferring certain liquid volumes and rinsing laboratory equipment

P.I.P. Integral Wash Bottle is made from polyethylene with a capacity of 250ml. It has a flexible body: squeeze it and after releasing, it will go back to the normal shape. Its ribbed screw-top cap facilitates opening and closing.

Integrated side spout eliminates possibility of contamination, dripping and leakage. P.I.P. Integral Wash Bottles are resistant to most of the common laboratory acids and chemicals.

GS1 Code	Capacity	Dimensions (mm)
6260807500514	250 ml	71 × 45 × 160



1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.