



Pole Ideal Pars Co.
www.medpip.com

Mega-Mix Rack

P.I.P. Mega-Mix Racks are available in two models of 150-well (for tubes with maximum diameter of 16mm) and 216-well (for tubes with maximum diameter of 12mm). Each rack has three movable plates, so that by interchanging the plates of the two models, users can place different sizes of test tubes in different heights into one rack. These racks are

made from sturdy polypropylene in different colors. Their design allows user to disassemble the rack to save space when necessary. Mega-Mix Racks are autoclavable and resistant to most of the common laboratory acids and chemicals, available in different colors for easier classification.

- Alphanumeric grid
- Robust structure
- Height - adjustable plates

GS1 Code	Model	Well Diameter	Dimensions (mm)	Wells/Plate	Wells/Rack
6260807501245	216-well	φ 12 mm	335 × 233 × 90	72 wells (12 × 6)	216 wells (18 × 12)
6260807501252	150-well	φ 16 mm	335 × 233 × 90	50 wells (10 × 5)	150 wells (15 × 10)



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran
Postal code:1576635714 P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159
Tel: 88765847 - 88765915 - 88545922-9 info@medpip.com www.medpip.com
@poleideal @poleideal +98 912 334 0197 Pole Ideal Pars Co.