




Pole Ideal Pars Co.
www.medpip.com

Microtube Combo Rack

Storing and organizing microtubes while keeping them cold

P.I.P. Microtube Combo Rack is designed in two models, suitable for accommodating different microtubes in one box which allows user to keep the samples cold. To do so, the box's lower container can be filled with crushed ice or water. This rack is made from

polypropylene which is autoclavable and freezable and resistant to most of the common laboratory acids and chemicals. Its microtube holder has alphanumeric grid for easy retrieval of the microtubes. Microtube Combo Rack has a transparent hinged lid with a latch.



GS1 Code	Model	Appropriate Tube (ml)	Number of Wells	Outer Dimensions (mm)	Qty./pack
6260807502761	69-well	0.2	30	176 × 95 × 55	5
		0.5	24		
		1.5 & 2	15		
6260807502778	54-well	0.2	12	176 × 95 × 55	5
		0.5	12		
		1.5 & 2	30		

1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

Postal code:1576635714

P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159

Tel: 88765847 - 88765915 - 88545922-9

info@medpip.com

www.medpip.com

@poleideal

@poleideal

+98912 334 0197

Pole Ideal Pars Co.