



Pole Ideal Pars Co.  
www.medpip.com

## Microtube Storage Rack

### Storing and organizing microtubes

P.I.P. Microtube Storage Racks are made from polypropylene, are autoclavable, freezable, and resistant to most of the common laboratory acids and chemicals. Microtube Storage Rack has alphanumeric grid for an easy retrieval of microtubes. These racks have

transparent hinged lids with latch, and their rectangular structure as well as the raised edges on the lid makes them safely stackable. They are available in 3 different models and various colors. The variety of colors facilitate classification.

GS1 Code	Model	Appropriate Microtube	Dimensions (mm)	Qty./pack
6260807500866	153 well (17 × 9)	0.2 ml	176 × 95 × 55	5
6260807500873	98 well (14 × 7)	0.5 ml	176 × 95 × 55	5
6260807500880	72 well (6 × 12)	1.5 ml	176 × 95 × 55	5



1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

Postal code:1576635714

P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159

Tel: 88765847 - 88765915 - 88545922-9

info@medpip.com

www.medpip.com

@poleideal

@poleideal

+98912 334 0197

Pole Ideal Pars Co.