



## MultiRack

### Storing and organizing different test tubes and microtubes

P.I.P. MultiRacks are available in 12 models, useful for storing and managing different kinds of tubes and microtubes in different laboratory departments. They are made from POM, submersible and resistant to most of the common laboratory acids and chemicals. The

handles of these racks facilitate carrying the samples. They are compact and lightweight, so they are ideal for placing inside incubators, refrigerators, freezers, under lab hoods and on bench tops. MultiRacks are provided in different colors.

GS1 Code	Model	No. of wells	Appropriate Tubes	Dimensions (mm)
6260807501924	L12 - 84	84	Tubes with diameters up to 13mm 12 × 100 mm tubes, round bottom tubes	293 × 115 × 65
6260807501894	S12 - 42	42	Tubes with diameters up to 13mm 12 × 100 mm tubes, round bottom tubes	175 × 115 × 65
6260807501900	L16 - 60	60	Tubes with diameters up to 16mm 16 × 100 mm tubes, 15ml conical tubes	293 × 115 × 65
6260807501870	S16 - 30	30	Tubes with diameters up to 16mm 16 × 100 mm tubes, 15ml conical tubes	175 × 115 × 65
6260807501863	L30 - 18	18	Tubes with diameters 25 to 30 mm 50 ml conical tubes	293 × 115 × 65
6260807501931	S30 - 9	9	Tubes with diameters 25 to 30 mm 50 ml conical tubes	175 × 115 × 65
6260807501955	LM - 60	60	1.5 and 2 ml microtubes	293 × 115 × 37
6260807501948	SM - 30	30	1.5 and 2 ml microtubes	175 × 115 × 37
6260807502211	LC - 84	84	Tubes with diameters up to 11 mm Round bottom tubes, cryotubes	293 × 115 × 37
6260807502228	SC - 42	42	Tubes with diameters up to 11 mm Round bottom tubes, cryotubes	175 × 115 × 37
6260807502266	LF - 39	39	15 and 50 ml conical tubes	293 × 115 × 65
6260807502235	LT - 72	72	Test tubes, round bottom tubes 12 × 100 mm and 16 × 100 mm tubes	293 × 115 × 65