



Outer Dimensions

Length : 440 mm
Width : 280 mm
Height with handle : 280 mm
Height without handle : 160 mm

Phlebotomy Tray

Safety transportation of samples and sampling tools

P.I.P. Phlebotomy Tray is specially designed to ensure that the sampling procedure is carried out quickly and specimens are handled safely and properly. It features a molded-in handle for easy transport and its unique bottom structure

allows the user to put the tray comfortably both on flat surfaces and most of the bed rails. P.I.P. Phlebotomy Tray is made from polypropylene, available in different colors which facilitate classification.

- Without lid, with deep compartments, preventing spillage and falling of equipment
- Customizing compartments using kits and dividers
- Special design for placing on horizontal poles (bed rails, trolleys, etc.)
- Lightweight and robust

GS1 Code	Model	Well diameter	Applications
6260807502082	General Kit	-	Transporting phlebotomy and sampling instruments, accommodating tube racks
6260807502099	Tube Kit	ϕ 44 mm (4 wells) ϕ 40 mm (4 wells) ϕ 21 mm (4 wells) ϕ 17 mm (14 wells) ϕ 13 mm (18 wells)	Holding different test tubes and sampling containers, suitable for transporting phlebotomy and sampling instruments, accommodating tube racks
6260807502105	Container Kit	ϕ 44mm (10 wells) ϕ 48mm (4 wells)	Holding different sample and stool containers, suitable for transporting phlebotomy and sampling instruments, accommodating tube racks

