



Pole Ideal Pars Co.
www.medpip.com

Plastic Beaker with Handle

Laboratory general use and handling liquids

P.I.P. Graduated Beakers are autoclavable, freezable, and made from transparent polypropylene. A proper and secure handle is embedded on this container to easily carry

and pour fluids. They are also resistant to most of the common laboratory acids and chemicals and are available in four different capacities.

GS1 Code	Capacity	Increments	Diameter × Height (mm)
6260807500217	250 ml	10 ml	φ 72 × 100
6260807500224	500 ml	10 ml	φ 85 × 124
6260807500231	1000 ml	10 ml	φ 96 × 180
6260807500248	2000 ml	20 ml	φ 126 × 210



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

Postal code:1576635714

P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159

Tel: 88765847 - 88765915 - 88545922-9

info@medpip.com

www.medpip.com

@poleideal

@poleideal

+98912 334 0197

Pole Ideal Pars Co.