



Pole Ideal Pars Co.
www.medpip.com




Plastic Kidney Dish

Receiving general medical and laboratory wastes

P.I.P. Plastic Graduated Kidney Dish is used in hospital general and emergency departments, operating rooms, dental clinics, etc. to receive soiled dressings and other medical waste before disposal. Made from polypropylene, these containers are autoclavable and resistant to most of the common laboratory acids and chemicals. Its capacity is 500cc and it has cc (every 100cc) and fluid ounce (every 4fl.oz.) graduations.

GS1 Code	Capacity (cc)	Capacity (fl.oz.)	Dimensions (mm)	Qty./Pack
6260807500767	500	16	220 × 107 × 46	150



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.