



## Pop-Up Rack

Carrying and organizing laboratory sample containers

P.I.P. Pop-Up Racks are made from polypropylene, are autoclavable, freezable, and resistant to most of the common laboratory acids and chemicals. In addition to stabilizing sample containers and tubes, these

racks provide easy storage and organization of containers and test tubes. They are available in different colors. This product is foldable to save space while stacking and autoclaving.

GS1 Code	Model	Dimensions (mm)	Number of Wells	Well Diameter
6260807501085	Rack 58	325 × 242 × 65	12 (3 × 4)	φ 58 mm
6260807501092	Rack 48	325 × 242 × 65	20 (4 × 5)	φ 48 mm
6260807501108	Rack 43	325 × 242 × 65	20 (4 × 5)	φ 43 mm
6260807501115	Falcon Rack *	325 × 242 × 65	40 (5 × 8) + 28 (4 × 7)	φ 17 mm & φ 29 mm

\* With alphanumeric grid

