



## Rotamix PIT180

### Mixing samples with different viscosities

Rotamix PIT180 is one of the basic devices for all laboratories to mix liquid samples with different viscosities. The rotation speed is adjustable from 10 to 80 rpm. There are various rotisserie rotors with different clips for this tube rotator to grip test tubes from 1.5 to 50 ml. Rotamix PIT180 is

ideal for general applications in medical, research, and industrial laboratories due to supporting different tubes and adjustable rotation time and speed. Rotamix PIT180 can be used in cold environments or inside incubators (5-40°C).

- Ergonomic and beautiful design
- Noiseless rotation system
- Adjustable rotation speed and time
- Emergency stop (by pressing STOP key twice)
- LCD display safety in case of liquids and solvents spillage

Technical Data	Rotation speed	10-80 rpm
	Voltage	210-230 V
	Frequency	50 Hz
	Max. Power	40 W
	Display	LCD
	Timer range	99 : 59' : 59"
	Permissible ambient temperature	5 - 40°C
	Max. relative humidity	80%
	Weight (approx.)	8 Kg

GSI Code	Cat. No.	Dimensions
6265738000447	PIT180	520 mm × 225 mm × 250 mm

