

# Rotary Micropipette Stand (6 - Clamp)



Pole Ideal Pars Co.  
[www.medpip.com](http://www.medpip.com)

# Rotary Micropipette Stand (6 - Clamp)

- Made from ABS
- Highly stable
- Occupation small space
- Holding and organizing 6 micropipettes
- Rotary axis for facilitates retrieving micropipettes
- Appropriate for most of pipettes (different size and brand )
- Resistant to most of the common laboratory acids and chemicals



GS1 Code	Model	Base Diameter	Height (mm)
6260807500668	6-clamp	ø154	340

