

P.I.P. Round bottom tube is disposable which is available non-sterile in two models of transparent (polystyrene) and opaque (polypropylene). They are autoclavable and can be sterilized by Gamma ray and Ethylene Oxide. These tubes are greatly applicable in Immunoassay tests (RIA, EIA, etc) which can be used for storage and dilution. These tubes are designed in dimensions that can be directly put in many analyzers so that changing the tube will not be needed and prevents sample from spilling and being wasted while transferring.



1st floor , No.12, Naghdi St., Jahantab St.,  
Mottahari Ave. Tehran, Iran  
Tel: +98 21 88545922-9  
Fax: +98 21 88765561 & +98 21 88767159  
Postal code: 1576635714 P.O. Box: 15875-9483



FOR SAFETY BOX

[www.medpip.com](http://www.medpip.com)  
[info@medpip.com](mailto:info@medpip.com)

Cat. No.: 92146



**Pole Ideal Pars**

[www.medpip.com](http://www.medpip.com)

**Round bottom tube**



Pole Ideal Pars

## Round bottom tube

### Technical Information

In 2 models

- Opaque (polypropylene) ○
- Transparent (polystyrene) ○

- Disposable, non-sterile
- Storing laboratory samples
- Suitable for dilution and immunoassay tests

Model	Transparency	Height x Diameter *
PS	Transparent	75 x Ø 12
PP	Opaque	75 x Ø 12