



Pole Ideal Pars Co.
www.medpip.com

Simple Test Tube

Storage and transport of laboratory samples

These disposable tubes are made from transparent polypropylene and can be used for storing or transporting different samples. They are resistant to most of the common laboratory acids and chemicals. Simple test tubes are available in two sizes, with cap (push cap made from polypropylene and snap cap made from polyethylene) or without cap.



GS1 Code	Capacity	Material - Cap type	Dimensions (diameter x height)	Qty./pack
6260807500422	7 ml	-	φ 12 mm x 100 mm	200
6260807500439	9 ml	-	φ 16 mm x 100 mm	200
6260807500446	7 ml	PP - Push	φ 12 mm x 100 mm	200
6260807500453	9 ml	PP - Push PE - Snap	φ 16 mm x 100 mm	200



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.