



Pole Ideal Pars Co.
www.medpip.com


Slide Holder with Grooved Slot

Holding and organizing slides

P.I.P. Slide Holder is made from rigid polystyrene and is ideal for holding 21 slides. It can be used to carry slides in laboratories. User can easily

place slides inside this slide holder and the color-coded slides will be easily separated and completely visible in an orderly manner.

GS1 Code	Capacity	Dimensions (mm)	Qty./Pack
6260807500842	21 wells	220 × 61 × 17	6



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran
Postal code:1576635714 P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159
Tel: 88765847 - 88765915 - 88545922-9 info@medpip.com www.medpip.com
@poleideal @poleideal +98912 334 0197 Pole Ideal Pars Co.