



## Slide Staining Set

Organizing staining solutions and increasing precision and efficiency in staining

P.I.P. Slide Staining Set is ideal for staining in histology, cytology, hematology, pathology, and microbiology. These sets consist of dishes with attached hinged lids, Slide Staining Holder with handle for 25 slides, and a Metal Holder. Slide Staining Dishes and the Slide Staining Holders are made from PET which is resistant to common solvents in the histology laboratories (such as toluene and xylene). The capacity of Slide Staining Dishes is 300ml and it is deep enough for the total immersion of

the slides. Besides, its lid is specially designed to be completely shut in order to prevent the reagents from probable evaporation. The Slide Staining Holder for 25 slides has a handle designed for convenient and controlled grip and transfer. The Metal Holder for organizing these dishes is made from stainless steel, it is washable and resistant to common laboratory acids and solvents. In addition, since all of these items are lightweight, user can handle and transfer them effortlessly.

GS1 Code	Model	Contents	Dimensions (mm)
6260807501603	Histology/ Cytology Set	12 Slide Staining Dishes + 2 Slide Staining Holders + 12-place Metal Holder	Metal Holder: 720 × 130 × 110 Staining dishes: 125 × 55 × 115 Staining holder: 92 × 30 × 92
6260807501610	Hematology Set	3 Slide Staining Dishes + 1 Slide Staining Holder + 3-place Metal Holder	Metal Holder: 230 × 130 × 110 Staining dishes: 125 × 55 × 115 Staining holder: 92 × 30 × 92

\* For more information on the chemical resistance of polystyrene (PS), refer to our website at [www.medpip.com](http://www.medpip.com).