



Pole Ideal Pars Co.
www.medpip.com

Specimen Container with Spoon

Sampling and storing the samples

P.I.P. Specimen Containers with Spoon are disposable, made from polypropylene, with a capacity of 15ml and depth of 50mm, suitable for storing and carrying different samples such as tissue, mucus, semen etc., with leakproof screw-top lid and attached spoon. In order

to facilitate storing, it is available in different colors which can help laboratory personnel to choose and assign certain colors to each sample. This container is autoclavable and resistant to most of the common laboratory acids and chemicals.

| GS1 Code | Capacity | Dimensions (diameter × height) | Total Height | Spoon height | Qty./pack |
|---------------|----------|-----------------------------------|--------------|--------------|-----------|
| 6260807501054 | 15 ml | φ 20 mm × 61 mm | 78 mm | 45 mm | 1000 |



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran
Postal code:1576635714 P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159
Tel: 88765847 - 88765915 - 88545922-9 info@medpip.com www.medpip.com
@poleideal @poleideal ☎ +98 912 334 0197 Pole Ideal Pars Co.