



Pole Ideal Pars Co.
www.medpip.com

Staining Tray with Slide Holder

Staining and holding slides

P.I.P. Staining Tray is ideal for staining different slides of histology, cytology, pathology etc. It has two separate sections for staining and holding slides. The staining chamber has a capacity of 21 slides at the same time and the holder has 21 grooves to hold the slides;

these parts can also be used separately. One of the advantages of this tray is its discharge pipe which provides the tidiness of the lab. It is resistant to most of the common laboratory acids and chemicals.

GS1 Code	Dimensions
6260807500545	340 mm × 225 mm × 35 mm



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

Postal code:1576635714

P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159

Tel: 88765847 - 88765915 - 88545922-9

info@medpip.com

www.medpip.com

📍 @poleideal

📷 @poleideal

☎ +98 912 334 0197

🏢 Pole Ideal Pars Co.