



Pole Ideal Pars Co.  
www.medpip.com


## Stainless Steel Tissue Capsule

Keeping and processing tissues

P.I.P. Tissue Capsule is autoclavable and made from stainless steel. There are holes embedded on its lid and base to facilitate the fluid flow. It has a snap cap with 3

clamps. Tissue Capsule is resistant to most of the common laboratory acids and chemicals and has secure snap cap.

GS1 Code	Dimensions (height × diameter)	Qty./pack
6260807500323	10 mm × φ 38 mm	10



1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.