



Stool Container with Spoon

Storing and transferring stool samples

The disposable P.I.P. Stool Container with an appropriate design is easy to use by patients and is available in two capacities of 15 and 30 ml with screw or push cap. The spoon is not embedded at the center of the lid so that it can be used

more easily without the possibility of hand touch. Moreover, the container's shape makes it proper for labeling. P.I.P. stool containers are made from polypropylene and polyethylene.



GS1 Code	Model	Spoon Type	Spoon Height	Dimensions (diameter × depth)	Capacity	Qty./pack
6260807500354	Small	Push	22 mm	25 mm × φ 32 mm	15 ml	2000
6260807500347	Large	Screw	32 mm	32 mm × φ 40 mm	30 ml	1200

