




Pole Ideal Pars Co.  
www.medpip.com

## Tak Rack Cryotube

Storing and organizing microtubes and cryotubes with 8 removable segments

P.I.P. Tak Rack Cryotube is designed with 8 segments to hold 64 x 1.5 or 2ml microtubes / cryotubes. Wells of this rack have a center depression to

stabilize conical tubes and a cryogenic vial lock for one-handed cap removal.



GS1 Code	Appropriate tube	Number of Wells	Outer Dimensions (mm)
6260807502709	1.5 & 2 ml (Cryo and microtubes)	64	140 x 140 x 55

1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran  
Postal code:1576635714 P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159  
Tel: 88765847 - 88765915 - 88545922-9 info@medpip.com www.medpip.com  
@poleideal @poleideal +98912 334 0197 Pole Ideal Pars Co.