



Pole Ideal Pars Co.
www.medpip.com


Tak Rack Microtube

Storing and organizing microtubes with 8 removable segments

P.I.P. Tak Rack Microtube has 8 reversible segments that 64 x 0.5ml microtubes by using one side of the rack or 96 x 0.2ml microtubes by using the other side can be placed in this

rack (they can also be combined). Another advantage of this product is the ease of use and it stops opened microtube lids from interfering with others.

GS1 Code	Appropriate tube	Number of Wells	Outer Dimensions (mm)
6260807502716	0.2	96	140 x 140 x 55
	0.5	64	



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

Postal code:1576635714

P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159

Tel: 88765847 - 88765915 - 88545922-9

info@medpip.com

www.medpip.com

@poleideal

@@poleideal

+98912 334 0197

Pole Ideal Pars Co.