



Pole Ideal Pars

lab dancer

LD8809

Test tube Shaker LD8809 model is a new product of P.I.P. with special features:

- With modern and efficient design
- Designed for mixing small amounts of
- Small, compact and reliable
- For small containers with maximum capacity of 30ml (e.g. test tubes, centrifuge tubes, and microtubes)
- With constant speed of 2600rpm
- Highest mixing Quality
- With a lightweight power supply and 12w power



lab dancer

Technical Specification:

Technical characteristic description	Amount
Type of motion	Orbital
Orbital motion amplitude	4.5 mm
Solution weight limit (including the container)	100 g
Device weight	550 g
Dimensions (height x diameter)	100mm x 70mm
Environment temperature limit	5-40 °C
Approximate humidity of environment	80%
Voltage	100-240 V
Frequency	50-60 Hz
Constant speed	2600 rpm

