



Pole Ideal Pars Co.
www.medpip.com

Applicator Stick (Without Cotton)

Mixing or blending substances

P.I.P. Applicators are used in biological sampling or for mixing and blending samples. They are available in two PS and PP models which are cotton free. Their proper length provides reachable and easy sampling. PS model is

made from polystyrene which is rigid and brittle. PP model is made from polypropylene which is soft and flexible. It is resistant to most of the common laboratory acids and chemicals.



GS1 Code	Model	Weight	Dimensions (length × diameter)	Qty./Pack
6260807500552	Brittle- PS	0.7 g	150 mm × ϕ 2.5 mm	500
6260807500569	Flexible - PP	0.7 g	150 mm × ϕ 2.5 mm	500



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.