


Chain Rack

Storing and organizing test tubes

P.I.P. Chain Rack has a unique linkable design that enables the possibility of multiple configurations to form different shaped tube racks. It holds 5 ml tubes on one end or 15 ml tubes (15-17 mm of diameter) on the

other. User can mix and match the well sizes to create one rack that holds multiple sized tubes in a variety of favorite and functional configurations.

- Saving bench space with the flexibility of holding multiple tube sizes in one configurable rack
- Creating larger tube rack configurations by linking together multiple Chain Racks
- Made from polypropylene



GS1 Code	Appropriate Tube	Outer Dimensions	Qty./pack
6260807502693	5 & 15 ml φ 15-17 mm	97 mm × φ 27 mm	42

Different Layouts:



Phlebotomy Tray (page 72)



Biojar (page 66)



Sample Collection and Transportation Case (page 73)



Cool Box (page 70)



Injection and Dressing Tray (page 71)

