



Pole Ideal Tajhiz Co.
www.medpit.com

Rocker Shaker

Model PIT16D



Rocker Shaker

Rocker Shaker PIT16D is designed for general and specialized applications. This device is a basic equipment in all medical, petroleum, food analysis, educational, and chemical laboratories. Due to appropriate capacity, this shaker makes it possible to blend and combine a wide range of liquids inside different containers which meets the requirements of many related tests



Features

- Digital display
- Balanced while running
- Noiseless mechanical system
- Ergonomic and beautiful design
- Rotation speed from 5 to 80 rpm
- Possibility of connecting to computer
- Blending liquids with rocking movement
- Safety of the sensitive parts against possible spillage of the fluids being tested

Rocker Shaker PIT16D Technical Data

Motion Type	rocking
Motion range of tray (Rocking angle fixed)	16°
Max. symmetrical loading capacity (kg)	5.0
Max. running time	Unlimited
Function type	Time limited or unlimited time
Speed range (rpm)	5-80
Speed display	Digital
Timer	0 - 99:59'/∞
Timer display	Digital
Outer dimensions (mm)	420 x 360 x 130
Max Input / Output power (W)	10-45
Voltage (VAC)	210-230
Frequency (HZ)	50
Weight (approx.) (kg)	8.8
Ambient permissible temperature (C°)	5-50



وزارت صنعت، معدن و تجارت
پولایدهال پارس



ACCREDITED



CERTIFIED
ISO 13485:2004



CERTIFIED
ISO 9001:2015

Unit11, No.13, Hojjatdoust st., South Felestin St., Keshavarz Blvd., Tehran, Iran

Tel:+98 21 88994482-3

Fax:+98 21 88994484

Postal Code:1416694334

Sales Office: +98 21 88545922

info@medpit.com / www.medpit.com

@poleideal

09123340197

@poleideal

Pole Ideal Pars Co.