



Pole Ideal Pars Co.
www.medpip.com

Daily – Weekly Pill Reminder

Storing and carrying pills during a week

P.I.P. Daily – Weekly Pill Reminder is a simple, lightweight, and easily applicable product by which the user can organize the pills and take them on time, since it has 7 individual sections for seven days of the week, each of which has 4


compartments. This Pill Reminder is suitable for those who need to take a number of pills and capsules during the day or are apt to forget and also for carers.

- Weekly classification of pills using the 7 compartments
- Compartments with separate lids for morning, noon, evening, bedtime
- Bottoms to facilitate the removal of each seven compartments from the main body
- Separate compartments usable for travel, work, etc.
- See-through compartments
- Markings on the main body and lids (7 days of the week, 4 times of the day)
- Marked with braille



Note: Some medications require special storage conditions and some cannot be kept in the same container as others.

GS1 Code	Model	Dimensions (mm)
6260807502808	daily-weekly	210 × 125 × 25



1st Floor, NO.12, Naghdi St., Jahantab St., Motahari Ave., Tehran, Iran

Postal Code: 1576635714

PO.Box: 15875 - 9483

✉ Info@medpip.com

🏢 Pole Ideal Pars Co.

☎ +98 9123340197

✈ @poleideal

📷 poleideal

Tel: +98 21 88545922-9

Fax: +98 21 88765561

+98 21 88767159