



Pole Ideal Pars Co.  
www.medpip.com

## Pipette and Thermometer Stand

Holding different glass pipettes and thermometers

P.I.P. Pipette and Thermometer Stand is steady, autoclavable, made from polypropylene, and resistant to most of the common laboratory acids and chemicals. It is available in two models of 25 and 50-well. The 25-well model has wells with 8mm diameter, which is appropriate for 0.1 to 3 ml glass pipettes and

1 to 8 ml volumetric pipettes. Besides, all glass thermometers that have diameters less than 8mm can be placed in this model. The 50-well model has wells with 16mm diameter, which is appropriate for 5 to 25ml glass pipettes and 7 to 50ml volumetric pipettes. Both models can be collapsed and reconnected for saving space.

GS1 Code	Model	Well Diameter	Appropriate Pipette	Appropriate Volumetric Pipette	Dimensions (mm)
6260807501269	25-well	φ 8 mm	0.1 - 3 ml	1 - 8 ml	145 × 137 × 255
6260807501276	50-well	φ 16 mm	5 - 25 ml	7 - 50 ml	210 × 137 × 255



1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran  
Postal code:1576635714 P.O. Box:15875-9483 FAX:+98 21 88765561 /+98 21 88767159  
Tel: 88765847 - 88765915 - 88545922-9 info@medpip.com www.medpip.com  
@poleideal @poleideal ☎ +98 912 334 0197 Pole Ideal Pars Co.