



Pole Ideal Pars Co.
www.medpip.com

Plastic Tongue Depressor

Examination of mouth and throat, blending drugs and ointments

P.I.P. disposable Tongue Depressor is made from polystyrene in two models of white and transparent with smooth surface and rounded edges and no unpleasant odor or taste. What makes this product special is its curved body which facilitates looking inside patient's

mouth, since physician's hand will be in a lower level than patient's mouth. More rigidity and flexibility and not having a mucosa-damaging rough surface are some other advantages of this tongue depressor. It can also be used to apply ointment, mixing drugs etc.



GS1 Code	Capacity	Dimensions (length x width)	Qty./pack
6260807500798	6 g	150 mm x 18 mm	100
6260807501344	6 g	150 mm x 18 mm	150

* For more information on the chemical resistance of polystyrene (PS), refer to our website at www.medpip.com.



1st floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.