



Sample Collection Funnel

Ease of collecting saliva, sputum, sperm and...



P.I.P sample collection funnel is ideal for collecting saliva, sputum, sperm, etc. This product can be connected to different types of tubes. It is autoclavable. P.I.P. funnel allows user to pour the sample into the tube easily, without spreading & contamination. By eliminating the steps of pipetting and transferring sample from containers to tubes, it allows easy processing of the sample.

- Ease of sample collection and processing
- Curved edge
- Made of polypropylene (PP)
- Suitable for collecting samples for pharmacogenomic tests, PCR, SNP genotyping, next generation sequencing, etc.

GS1 Code	Model	Size	Dimension (mm)
6260807503843	With 5ml conical tube	Small	28x45x90
6260807503775	With 5ml conical tube (Amber)	Small	28x45x90
6260807503836	With 15ml conical tube	Small	28x45x144
6260807503805	With 15ml conical tube (Amber)	Small	28x45x144
6260807503799	With 25ml conical tube	Larg	29x53x120
6260807503812	With 25ml conical tube (Amber)	Larg	29x53x120
6260807503768	With 50ml conical tube	Larg	29x53x160
6260807503829	With 50ml conical tube (Amber)	Larg	29x53x160
6260807503782	With Sampling tube with spoon	Larg	29x53x96
6260807503850	with test tube 16 ×100 mm	Small	28x45x126
6260807503867	with test tube 16 ×54 mm	Small	28x45x78

